

Bagmatic SL



Barcode scanner



High capacity battery



Radio data transfer

Blood donation genius

- Most reliable and stable technology worldwide
- Complete control of all donations
- 3 ways of safe data transfer
- Analysis of donation data
- Entirely free programmable and adaptable to your needs
- Free configuration of the software

... and many other features!

Robust transport case



Chargematic



Technical data

Power Supply	12 V DC
Power Consumption	0,15 A/ 0,35 A
Fuse	T 0.8A ; 250 V
Battery charger	FRIWO FW 7118M/12
Battery charging time	6 hours
Surrounding temperature	10°C - 30°C
Humidity	maximum 85 %
Dimensions	365x190x135 mm
Weight	5,1 kg