

BT100F dispensing intelligent peristaltic pump

- Color touch screen, making the information more intuitive and the operation more convenient;
- Intelligent expert operating system can guarantee smooth operation even the operator never used such kind of machine before;
- Control technology of motor rotation angle precision, reducing the dispensing errors under traditional time controlling distribution pattern due to motor slow start;
- The expert operating system can display the tube specifications, materials, pump head types and other information, guiding the operator set peristaltic pump correctly according to the input flow, fluid volume, working hours and other data;
- Intelligent temperature control technology greatly reduces the operating noise of peristaltic pump;
- It has three kinds of distribution patterns, time dispense mode, fluid volume dispense mode, copying dispense mode. It can meet different needs;
- The function of accumulated distribution of liquid can record total liquid volume;
- Infrared remote control function, it can be non-contact control;
- BT100F dust-proof of stainless steel case design, suitable for general industrial environments;
- For pump head: YZ15,2*YZ15, 4*YZ15,YZ25,2*YZ25, DG6-1,DG6-2,DG6-4,DG6-8,DG6-12, DG10-1,DG10-2,DG10-4,DG10-8,DG10-12;
- Flow range:0.00016~575(ml / min).

