


	SEE BELOW	
WORKING RANGE (µL)	50 – 250 µL	
MATERIAL	PS	
	MM	INCH
A, BASE LENGTH	127,8	5,03
B, BASE WIDTH	85,5	3,35
C, TOP LENGTH	123,3	4,84
D, TOP WIDTH	81,2	3,20
E, FLANGE LONG SIDE	2,4	0,08
F, OVERALL HEIGHT	14,4	0,55
G1, DIA. TOP	Ø7	Ø0,28
G2, DIA. BOTTOM	Ø6,2	Ø0,24
H, A - I LOCATION (X)	14,4	0,55
J, A - I LOCATION (Y)	11,4	0,45
K, WELL TO WELL SPACING	9	0,35
L, WELL DEPTH	11,4	0,45
M, BOTTOM OF WELL DISTANCE	3,2	0,13
N, FLANGE SHORT SIDE	2,4	0,10

CATALOG NO	SURFACE	COLOR	STERILE	WITH LID*	UNITS/PK/CAS
436006	AMINO	CLEAR	-	-	5/30
436007	AMINO	WHITE	-	-	5/30
436008	AMINO	BLACK	-	-	5/30
436014	STREPTAVIDIN	CLEAR	-	-	1/15
436015	STREPTAVIDIN	WHITE	-	-	1/15
436016	STREPTAVIDIN	WHITE	-	-	1/15
436024	NICKEL-CHELATE	CLEAR	-	-	1/15
436026	NICKEL-CHELATE	BLACK	-	-	1/15
436027	NICKEL-CHELATE	WHITE	-	-	1/15
436032	GLUTAHIONE	CLEAR	-	-	1/15
436033	GLUTAHIONE	WHITE	-	-	1/15
436034	GLUTAHIONE	BLACK	-	-	1/15

NOTES:

- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
- THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATION PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



ThermoFisher
SCIENTIFIC

Roskilde Site · Kamstrupvej 90 · P.O.Box 280
DK-4000 Roskilde · Denmark
Tel +45 4631 2000 · Fax +45 4631 2042
E-mail: info.nunc@thermofisher.com

www.nuncbrand.com

SUBJECT: F96 MICROWELL PLATE IMMOBILIZER, PS		DRAWN BY/DATE:	CHANGE ORDER NO.:	VALID FROM:	APPROVED BY/DATE:	FORMAT:
ITEM NUMBER:	MATERIAL:	DRAWING NO.:	VERSION:	PAGE TYPE:	SCALE:	METRIC 1:1 1/8"=1"

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC