

## WG600S variable speed peristaltic pump

- stainless steel case
- buttons control the speed, rotating, start/stop, full-speed, easy to operate;
- Suitable for pump head : YZ35, 2×YZ35
- Flow range : 800~13000 (ml/min)



### Technical parameters

- Speed range:60~600 rpm
- Speed resolution: 1 rpm
- Supply Voltage: AC220V±10% 50/60HZ
- Hz V: less than 300W
- Working Environment Temperature: 0 ~ 40 °C
- Relative humidity: <80%
- Dimensions: (L × W × H): 340mm×235mm×175mm
- Weight: 13kg

External control functions:

external control start and stop , external control orientation, Isolated, TTL electrical level input optional; 0-5V/10V, 4-20m optional;